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The fundafom® destroys the foam where it forms: inside the reactor itself. It consists of a drive unit and a separator disc stack, flanged onto the top cover of a closed reactor.

Depending on size, the disc stack carries several conical discs with radial guides welded to their lower surface. The motor rotates the stack. Rising foam entering the stack is broken up or burst by centrifugal force. The released gas is discharged free of droplets through the hollow shaft and the gas outlet nozzle, while the liquid phase is thrown back into the reactor.

This largely eliminates the need for chemical anti-foam additives. In sterile applications this matters: without additives the surface tension of the medium stays lower, which means a higher oxygen transfer rate and better growth of bacteria, yeasts, cells and mycelia.

In many cases the working volume of the reactor can be increased considerably. Existing plants have achieved production increases of between 10 and 80 %.

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Sizes	• 00 to 8
Materials	• Product contact: AISI 316L stainless steel • Special materials such as Hastelloy, Monel, titanium on request
Sealing materials	• EPDM (standard) • PTFE or FKM for special applications
Shaft seal	• Double-acting mechanical seal • Sterile operation: double axial mechanical seal with a sterile barrier medium under overpressure • Non-sterile operation: barrier chamber, level monitoring on request
Bearings	• Rolling bearings, vibration-free and quiet running
Drive	• Directly driven three-phase motor • Explosion-proof design available • Foam sensing probe for intermittent foam available
Mounting	• On the vertical centreline of the reactor cover • Off-centre where the centre is occupied • Inclination 0 to 15° from vertical
Media	• Foaming gas/liquid mixtures
Operating mode	• Pressure or vacuum
Air and gas flow rate	• 0.5 to 100 m³/min depending on size
Drive rating	• 0.25 to 90 kW depending on size • Guide values; the exact figure follows from the type of foam and the operating viscosity
Sizing	• Based on the volume of gas in operating cubic metres that has to leave the reactor per minute

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Guide values. The drive rating depends on the type of foam and the operating viscosity and is calculated once these are known.

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Fermentation

Chemical anti-foam additives often affect the growth of yeasts, cells, bacteria and mycelia and have to be separated from the product again at considerable cost. The fundafom® largely eliminates them.

- Citric acid, submerged process
- Penicillin, tetracycline, streptomycin, oxytetracycline
- Cheese moulds (Penicillium roqueforti, Penicillium camemberti)
- Baker's yeast and fodder yeast from molasses
- Yeast from sulphite liquor, whey and starch-bearing waste water
- Enzyme production for detergents
- Steroids
- Single cell protein from methanol
- Bacterial vaccines

Chemical and Process Industry

In many processes agitation, gassing or chemical reaction produces unwanted heavy foam.

- Monomer stripping in PVC production
- Manufacture of organic pigments
- Methionine
- Phosphoric acid for fertiliser production
- Vacuum distillation of organic solvents
- Vacuum evaporation and vacuum filtration
- Venting of buffer tanks
- Neutralisation of distillation residues
- Rolling oil venting in aluminium rolling mills
- Degassing of crude oil

Benefit

- Higher output through a larger working volume
- Recovery of the gases
- No loss of product
- No fouling or blocking of vacuum systems by foam carry-over

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The required size follows from the volume of air or gas in operating cubic metres that has to leave the reactor per minute.

For sizing we offer a flow simulation (CFD).

For that we need the completed questionnaire together with a drawing of the cover of the reactor concerned.

The questionnaire is available on request:
sales@liquitec.ch.

Mass- und Konstruktionsänderungen vorbehalten. Revisionsdatum: 09/26
Specifications and design are subject to change without notice. Date of revision: 09/26